
CP-31398 dihydrochloride

产品编号: D56150

CAS: 1217195-61-3

分子式: C₂₂H₂₈Cl₂N₄O

纯度: 99.53%

InChi Key: InChIKey=WWIFDJIOGCGSBS-IVKCLRODSA-N

Smiles: CN(C)CCCNC1=C2C=CC=CC2=NC(/C=C/C3=CC=C(OC)C=C3)=N1.[H]Cl.[H]Cl

外观: 固体粉末

作用通路: Apoptosis

作用靶点: MDM-2/p53; p53

溶解性: DMSO: 6 mg/ml (13.78 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: CP 31398 dihydrochloride is a p53 stabilizer which stabilizes the active conformation of p53 and promotes p53 activity in cancer cell lines with mutant or wild-type p53.