
GS-9667

产品编号: D56224

CAS: 618380-90-8

分子式: C₂₁H₂₄FN₅O₄S

纯度: 99.76%

InChi Key: InChIKey=IZRXENCTXNMAMI-DIJFLQFKSA-N

Smiles: O[C@H]1[C@H](N2C=3C(N=C2)=C(N[C@H]4[C@H](O)CCC4)N=CN3)O[C@H](CSC5=C(F)C=CC=C5)[C@H]1O

外观: 固体粉末

作用通路: GPCR/G Protein;Neuroscience

作用靶点: Adenosine Receptor

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: GS-9667, a selective and partial agonist of the A(1) adenosine receptor (AR), represents an effective therapy for Type 2 diabetes (T2DM) and dyslipidemia via lowering of free fatty acids (FFA).