
O4I1

产品编号: D56260

CAS: 175135-47-4

分子式: C₁₆H₁₅NO₂

纯度: 98.33%

InChi Key: InChIKey=QJEJLARVLKHYFW-UHFFFAOYSA-N

Smiles: COc1ccc(cc1)COc1ccc(cc1)CC#N

外观: 固体粉末

作用通路: Membrane transporter/Ion channel

作用靶点: OCT

溶解性: DMSO:48 mg/mL (189.5 mM);Ethanol:1 mg/mL (3.95 mM);Water:<1 mg/mL

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: O4I1 is an effective Oct3/4 inducer.