
SM-324405

产品编号: D56265

CAS: 677773-91-0

分子式: C₁₉H₂₃N₅O₄

纯度: 96.44%

InChi Key: InChIKey=MXILFZWMVNDIOZ-UHFFFAOYSA-N

Smiles: CCCCOc1nc(N)c2[nH]c(=O)n(Cc3cccc(CC(=O)OC)c3)c2n1

外观: 固体粉末

作用通路: Immunology/Inflammation

作用靶点: TLR

溶解性: DMSO: 9 mg/mL (23.35 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: SM 324405 is an agonist of TLR7.