

---

## Acetyltrialanine

产品编号: D56353

CAS: 19245-85-3

分子式: C<sub>11</sub>H<sub>19</sub>N<sub>3</sub>O<sub>5</sub>

纯度: 99.53%

InChi Key: InChIKey=DRYOODAJROGPQO-ACZMJKKPSA-N

Smiles: [C@H](C(N[C@H](C(O)=O)C)=O)(NC([C@@H](NC(C)=O)C)=O)C

外观: 固体粉末

作用通路: Others

作用靶点: Others

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Acetyltrialanine is a dipeptide compound that binds at two sites on the Tb<sup>3+</sup>-pancreatic elastase complex and can be used as a nitrogen source.