
IRAK-4 protein kinase inhibitor 2

产品编号: D56480

CAS: 301675-24-1

分子式: C₁₄H₁₀N₄O₃

纯度: 98%

InChi Key: InChIKey=QPKROMNIEAEFQA-UHFFFAOYSA-N

Smiles: [O-][N+](=O)C1=CC(=CC=C1)C(=O)NC1=NC2=C(N1)C=CC=C2

外观: 固体粉末

作用通路: Immunology/Inflammation; NF-κB

作用靶点: IRAK

溶解性: DMSO : 10 mg/mL (35.43 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: IRAK-4 protein kinase inhibitor 2 is a potent interleukin-1 receptor-associated kinase-4 (IRAK-4) inhibitor(IC₅₀ = 4 μM).