
Thrombin inhibitor 5

产品编号: D56602

CAS: 328108-09-4

分子式: C₁₁H₉FN₄O₃

纯度: 99.4%

InChi Key: InChIKey=CYAQTWKBNHAHHY-UHFFFAOYSA-N

Smiles: COC(=O)c1nnc(NC(=O)c2ccc(F)cc2)[nH]1

外观: 固体粉末

作用通路: Others

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Thrombin inhibitor 5 (compound 385) is a thrombin inhibitor with IC₅₀s ranging from 0.1 μM to 1 μM. Thrombin inhibitor 5 can be used for venous thromboembolism research.