
Rentiapril racemate

产品编号: D56611

CAS: 72679-47-1

分子式: C₁₃H₁₅NO₄S₂

纯度: 98%

InChi Key: InChIKey=BSHDUMDXSRLRBI-UHFFFAOYSA-N

Smiles: O=C(O)C1N(C(=O)CCS)C(SC1)C=2C=CC=CC2O

外观: 固体粉末

作用通路: Metabolism

作用靶点: Angiotensin-converting Enzyme (ACE)

溶解性: DMSO: 50 mg/mL (159.55 mM, Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Rentiapril racemate (SA-446 racemate) is the racemate of Rentiapril. Rentiapril racemate has anti-inflammatory activity and may be used in glaucoma research.