
CTX-0124143

产品编号: D56615

CAS: 423731-64-0

分子式: C₁₇H₁₃FN₂O₃S

纯度: 99.16%

InChi Key: InChIKey=CMGWGSMHPCWLOH-UHFFFAOYSA-N

Smiles: FC1=C(C=CC=C1)C(=O)NNS(=O)(=O)C1=CC2=CC=CC=C2C=C1

外观: 固体粉末

作用通路: Chromatin/Epigenetic

作用靶点: Histone Acetyltransferase

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: CTX-0124143 is a histone acetyltransferase inhibitor with an IC₅₀ value of 1.0 μM for KAT6A. CTX-0124143 can be used to study cellular senescence.