
rac-NBI-74330

产品编号: D56634

CAS: 473722-68-8

分子式: C₃₂H₂₇F₄N₅O₃

纯度: 97.85%

InChi Key: InChIKey=XMRGQUUDUVGRCBS-UHFFFAOYSA-N

Smiles: O=C1C=2C=CC=NC2N=C(N1C3=CC=C(OCC)C=C3)C(N(C(=O)CC4=CC=C(F)C(=C4)C(F)(F)F)CC=5C=NC=CC5)C

外观: 固体粉末

作用通路: Autophagy;GPCR/G Protein;Immunology/Inflammation

作用靶点: CXCR

溶解性: DMSO: 27.3mg/mL (45mM;)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: rac-NBI-74330 is an effective and selective CXCR3 antagonist.