
UVI 3003

产品编号: D56659

CAS: 847239-17-2

分子式: C₂₈H₃₆O₄

纯度: 99.94%

InChi Key: InChIKey=APJSHECCIRQQDV-ZRDIBKRKSA-N

Smiles: CCCCCOC1=C(C2=CC(/C=C/C(O)=O)=CC=C2O)C=C3C(C(C)(CCC3(C)C)C)=C1

外观: 固体粉末

作用通路: Autophagy; Metabolism

作用靶点: Autophagy; Retinoid Receptor

溶解性: DMSO:100 mg/mL (229.05 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: UVI 3003 is a highly selective antagonist of the retinoid X receptor. UVI 3003 inhibits Xenopus and human RXR α in Cos7 cells (IC₅₀s: 0.22 and 0.24 μ M, respectively).