
Enasidenib mesylate

产品编号: D56698

CAS: 1650550-25-6

分子式: C₂₀H₂₁F₆N₇O₄S

纯度: 98%

InChi Key: InChIKey=ORZHZQZYWXEDDL-UHFFFAOYSA-N

Smiles: O=S(=O)(O)C.FC(F)(F)C1=NC(=CC=C1)C=2N=C(N=C(N2)NCC(O)(C)C)NC=3C=CN=C(C3)C(F)(F)F

外观: 固体粉末

作用通路: Metabolism

作用靶点: Dehydrogenase

溶解性: DMSO: 80 mg/mL(140.48 mM,Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Enasidenib mesylate (AG-221 mesylate) is a potent and selective IDH2 mutase inhibitor that promotes differentiation of leukemia myeloid cells for the treatment of acute myeloid leukemia.