

---

## MK-8245 analog

产品编号: D56713

CAS: 1030612-87-3

分子式: C<sub>16</sub>H<sub>15</sub>BrFN<sub>7</sub>O<sub>3</sub>S

纯度: 99.71%

InChi Key: InChIKey=SVTBUTPWEOUGGT-UHFFFAOYSA-N

Smiles: OC(Cn1nnc(-c2nnc(N(CC3CCC3Oc(cc(cc3)F)c3Br)s2)n1)=O

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO: 55 mg/ml (113.57 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: MK-8245 analog is an analog of MK-8245 which is a liver-targeted Inhibitor of Stearoyl-CoA Desaturase (SCD). It is useful for the Treatment of Dyslipidemia and Diabetes.