

---

## Capzimin

产品编号: D56779

CAS: 2084868-04-0

分子式: C<sub>30</sub>H<sub>24</sub>N<sub>6</sub>O<sub>2</sub>S<sub>4</sub>

纯度: 98.31%

InChi Key: InChIKey=RNEOHKZPZKEZCQ-UHFFFAOYSA-N

Smiles: O=C(NCCCC1=NC=CS1)C2=CN=C3C(SSC4=CC=CC5=CC(=CN=C45)C(=O)NCCC6=NC=CS6)=CC=C3=C2

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO: 50.1 mg/mL (79.7 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Capzimin is a selective inhibitor of proteasome isopeptidase Rpn11.