
WIKI4

产品编号: D56804

CAS: 838818-26-1

分子式: C₂₉H₂₃N₅O₃S

纯度: 99.32%

InChi Key: InChIKey=RNUXIZKXJOGYQP-UHFFFAOYSA-N

Smiles: COC1=CC=C(C=C1)N1C(SCCCN2C(=O)C3=CC=CC4=C3C(=CC=C4)C2=O)=NN=C1C1=CC=NC=C1

外观: 固体粉末

作用通路: Chromatin/Epigenetic; DNA Damage/DNA Repair; Cytoskeletal Signaling; Stem Cells

作用靶点: PARP; Wnt/beta-catenin

溶解性: DMSO:7 mg/mL; Water:<1 mg/mL; Ethanol:<1 mg/mL

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: WIKI4 is a potent inhibitor of Wnt/ β -catenin signaling and tankyrase 2 (TNKS2).