
Fendosal

产品编号: D56810

CAS: 53597-27-6

分子式: C₂₅H₁₉NO₃

纯度: 96%

InChi Key: InChiKey=HAWWPSYXSLJRBO-UHFFFAOYSA-N

Smiles: OC(c(cc(cc1)-n2c(-c3ccccc3)cc-3c2CCc2c3cccc2)c1O)=O

外观: 固体粉末

作用通路: Immunology/Inflammation

作用靶点: Immunology/Inflammation related

溶解性: DMSO: 22.5 mg/mL (58.99 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Fendosal is a non-steroidal anti-inflammatory agent. It is 6.9 to 9.5 times more active than aspirin in the prophylactic and therapeutic adjuvant-induced polyarthritis models of chronic inflammation.