
Pipamperone

产品编号: D56812

CAS: 1893-33-0

分子式: C₂₁H₃₀N₃O₂

纯度: 99.94%

InChi Key: InChIKey=AXKPFOAXAHJUAG-UHFFFAOYSA-N

Smiles: NC(C1(CCN(CCCC(c(cc2)ccc2F)=O)CC1)N1CCCCC1)=O

外观: 固体粉末

作用通路: GPCR/G Protein; Neuroscience

作用靶点: 5-HT Receptor; Dopamine Receptor

溶解性: DMSO : 100 mg/mL (266.33 mM; Need ultrasonic); 0.1 M HCL : 25 mg/mL (66.58 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Pipamperone binds 5-HT_{2A} closely as receptor.