
LY404039

产品编号: D56825

CAS: 635318-11-5

分子式: C₇H₉NO₆S

纯度: 99.03%

InChi Key: InChiKey=AVDUGNCTZRCAHH-NFKAWTSISA-N

Smiles: C1[C@]([C@@H]2[C@H]([C@@H]2S1(=O)=O)C(=O)O)(C(=O)O)N

外观: 固体粉末

作用通路: Endocrinology/Hormones; Neuroscience

作用靶点: Glucocorticoid Receptor; GluR

溶解性: DMSO:10mM

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: LY404039 is an effective agonist of recombinant human mGlu2(Ki =149 nM)/mGlu3(Ki =92 nM).