
KW-2478

产品编号: D56841

CAS: 819812-04-9

分子式: C₃₀H₄₂N₂O₉

纯度: 98.68%

InChi Key: InChIKey=VFUXSYAXEKYYMB-UHFFFAOYSA-N

Smiles: CCc1c(O)cc(O)c(C(=O)c2cc(OC)c(OCCN3CCOCC3)cc2)c1CC(=O)N(CCOC)CCOC

外观: 固体粉末

作用通路: Cytoskeletal Signaling; Metabolism

作用靶点: HSP

溶解性: DMSO:106 mg/mL (184.5 mM);Ethanol:3 mg/mL (5.22 mM);Water:<1 mg/mL

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: KW-2478 is a non-ansamycin HSP90 inhibitor with IC₅₀ of 3.8 nM.