
Ossirene

产品编号: D56878

CAS: 106566-58-9

分子式: C₂H₄Cl₃O₂Te·H₄N

纯度: 98%

InChi Key: InChIKey=DXZFFLRJVDZCMT-UHFFFAOYSA-N

Smiles: N.Cl[TeH]1(Cl)(Cl)OCCO1

外观: 固体粉末

作用通路: Apoptosis; Proteases/Proteasome; Immunology/Inflammation

作用靶点: Caspase; IL Receptor; Interleukin

溶解性: DMSO:198.7 mM

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: AS101, a potent in vitro and in vivo immunomodulator, is a novel inhibitor of IL-1beta converting enzyme.