
3'-Deoxyguanosine

产品编号: D56903

CAS: 3608-58-0

分子式: C₁₀H₁₃N₅O₄

纯度: 98.85%

InChi Key: InChIKey=OROIAVZITJBGSM-OBXARNEKSA-N

Smiles: Nc1nc2c(ncn2[C@@H]2O[C@H](CO)C[C@H]2O)c(=O)[nH]1

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO:10mg/ml(37.4mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: 3'-Deoxyguanosine is a ligand that can be complexed with enzymes, such as purine nucleoside phosphorylase, and receptors.