

---

## Ozagrel hydrochloride

产品编号: D56905

CAS: 78712-43-3

分子式: C<sub>13</sub>H<sub>12</sub>N<sub>2</sub>O<sub>2</sub>·HCl

纯度: 99.78%

InChi Key: InChiKey=CWKFWBJJNPGAM-IPZCTEOASA-N

Smiles: c1cc(ccc1Cn1ccnc1)/C=C/C(=O)O.Cl

外观: 固体粉末

作用通路: GPCR/G Protein; Immunology/Inflammation; Proteases/Proteasome

作用靶点: Prostaglandin Receptor; Thrombin

溶解性: DMSO:49 mg/mL (185.1 mM);Ethanol:6 mg/mL (22.66 mM);Water:48 mg/mL (181.3 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Ozagrel HCl is a selective thromboxane A<sub>2</sub> synthetase inhibitor used for the improvement of postoperative cerebrovascular contraction and accompanying cerebral ischaemia.