
L-655708

产品编号: D56922

CAS: 130477-52-0

分子式: C₁₈H₁₉N₃O₄

纯度: 98%

InChi Key: InChIKey=YKYOQIXTECBVBB-AWEZLNQCLSA-N

Smiles: C(OCC)(=O)C1=C2N(C=3C(C(=O)N4[C@]2(CCC4)[H])=CC(OC)=CC3)C=N1

外观: 固体粉末

作用通路: Membrane transporter/Ion channel;Neuroscience

作用靶点: GABA Receptor

溶解性: DMSO:5 mg/mL (14.65 mM,Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: L-655708 is a selective and highly potent GABAA receptor reverse agonist with a Ki value of 0.45 nM.