
H2L5186303

产品编号: D56928

CAS: 139262-76-3

分子式: C₂₆H₂₀N₂O₈

纯度: 97.23%

InChi Key: InChIKey=HZFPOTBCYPWQSH-DZDAAMPGSA-N

Smiles: OC(=O)\C=C/C(=O)NC1=CC=C(OC2=CC(OC3=CC=C(NC(=O)\C=C/C(O)=O)C=C3)=CC=C2)C=C1

外观: 固体粉末

作用通路: GPCR/G Protein

作用靶点: LPA Receptor

溶解性: DMSO : 200 mg/mL (409.5 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: H2L5186303 is an LPA2 antagonist and can be used in studies about asthma treatment.