
ELISA包被液 (1×) ELISA Coating Solution(1×)

产品编号: C04-01001

保存条件: Store at 2-8°C for one year.

产品介绍: 本品为PH9.6的即用型工作液, 用于ELISA (酶联免疫吸附实验) 中酶标板的抗原或抗体的包被。使用时取适量包被缓冲液将抗原或抗体溶液稀释至相应包被的浓度, 在ELISA酶标板每个反应孔中加入100ul, 4°C包被过夜, 或37°C包被1h。

PRODUCT SPECIFIC PUBLICATIONS

[IF=6.6] Guibin Pang. et al. Improving combination cancer immunotherapy by manipulating dual immunomodulatory signals with enzyme-triggered, cell-penetrating peptide-mediated biomodulators. BIOMATER SCI-UK. 2023 Dec;; Other ; . 38167881

[IF=4.047] Liu, Yu-Ke. et al. Comparative analysis of the anticoagulant activities and immunogenicity of HSC70 and HSC70TKD of Haemaphysalis flava. PARASITE VECTOR. 2022 Dec;15(1):1-12 Other ; . 36335395