
Tolimidone

产品编号: D56938

CAS: 41964-07-2

分子式: C₁₁H₁₀N₂O₂

纯度: 100%

InChi Key: InChIKey=HJQILFPVRNHTIG-UHFFFAOYSA-N

Smiles: CC1=CC=CC(OC2=CNC(=O)N=C2)=C1

外观: 固体粉末

作用通路: Angiogenesis; Tyrosine Kinase/Adaptors

作用靶点: Src

溶解性: DMSO:20.2 mg/mL (100 mM);Ethanol:4 mg/mL (20 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: MLR-1023 is a chemical compound which inhibits acid secretion in animal models and also acts as a bronchodilator in histamine-challenged animals.