

---

## cim0216

产品编号: D56951

CAS: 1031496-06-6

分子式: C<sub>21</sub>H<sub>21</sub>N<sub>3</sub>O<sub>2</sub>

纯度: 99.6%

InChi Key: InChIKey=KSEXDSJYVSEVGF-UHFFFAOYSA-N

Smiles: Cc1cc(NC(=O)C(N2CCCC3CCCCC23)c2cccc2)no1

外观: 固体粉末

作用通路: Membrane transporter/Ion channel

作用靶点: TRP/TRPV Channel

溶解性: DMSO: 30 mg/mL (86.35 mM; Need ultrasonic and warming); Ethanol: 17 mg/mL (48.93 mM; Need ultrasonic and warming)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: CIM0216 is a synthetic TRPM3 activator whose potency and apparent affinity greatly exceeds that of the canonical TRPM3 agonist.