
COH-SR4

产品编号: D56963

CAS: 73439-19-7

分子式: C₁₃H₈Cl₄N₂O

纯度: 99.67%

InChi Key: InChIKey=VKVJIWVUYNTBEZ-UHFFFAOYSA-N

Smiles: ClC1=CC(NC(=O)NC2=CC(Cl)=CC(Cl)=C2)=CC(Cl)=C1

外观: 固体粉末

作用通路: PI3K/Akt/mTOR signaling; Chromatin/Epigenetic

作用靶点: AMPK

溶解性: DMSO: 4 mg/mL (11.43 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: SR4 is a uncoupler of mitochondrial oxidative phosphorylation. SR4 modulating amp-dependent kinase (ampk)-mammalian target of rapamycin (mtor) signaling, and inhibiting proliferation of hepg2 hepatocarcinoma cells