
SRX3177

产品编号: D56967

CAS: 2241237-51-2

分子式: C₃₁H₃₂N₆O₄S

纯度: 98%

外观: 固体粉末

作用通路: PI3K/Akt/mTOR signaling; Cell Cycle/Checkpoint; Chromatin/Epigenetic

作用靶点: PI3K; CDK; Epigenetic Reader Domain

溶解性: DMSO: 20 mg/mL (34.21 mM, Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: SRX3177 is a novel and potent inhibitor of BRD4, PI3K and CDK. SRX3177 has potential anticancer activity.