
Choline Fenofibrate

产品编号: D56971

CAS: 856676-23-8

分子式: C₂₂H₂₈ClNO₅

纯度: 99.95%

InChi Key: InChIKey=JWAZHODZSADEHB-UHFFFAOYSA-M

Smiles: CC(C)(C(=O)[O-])Oc1ccc(cc1)C(=O)c1ccc(cc1)Cl.C[N+](C)(C)CCO

外观: 固体粉末

作用通路: DNA Damage/DNA Repair; Metabolism; Immunology/Inflammation; Neuroscience

作用靶点: COX; P450; PPAR

溶解性: DMSO:6.2 mg/mL, Need ultrasonic and warming

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Choline Fenofibrate is a choline formulation of fenofibrate, a synthetic phenoxy-isobutyric acid derivate and prodrug with antihyperlipidemic activity.