
Ropsacitinib

产品编号: D56978

CAS: 2127109-84-4

分子式: C₂₀H₁₇N₉

纯度: 99.85%

InChi Key: InChIKey=XPLZTJWZDBFWDE-OYOVHJISSA-N

Smiles: C(C#N)[C@@]1(C[C@H](C#N)C1)N2C=C(C=N2)C=3C=4N(C=C(N3)C=5C=NN(C)C5)N=CC4

外观: 固体粉末

作用通路: Tyrosine Kinase/Adaptors

作用靶点: Tyrosine Kinases

溶解性: DMSO: 10mg/mL (26.08.2mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Ropsacitinib is a selective TYK2 inhibitor, which binds to TYK2 catalytically active JH1 domain (IC₅₀: 17 nM).