
Lobenzarit sodium

产品编号: D56988

CAS: 64808-48-6

分子式: C₁₄H₈ClNNa₂O₄

纯度: 98%

InChi Key: InChIKey=QPJBONAWFAURGB-UHFFFAOYSA-L

Smiles: [O-]C(c(cccc1)c1Nc(cc(cc1)Cl)c1C([O-])=O)=O.[Na+].[Na+]

外观: 固体粉末

作用通路: Others

作用靶点: Others

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Lobenzarit sodium is an agent with the activity of antirheumatic and antioxidative.