

---

## Lalistat 2

产品编号: D57032

CAS: 1234569-09-5

分子式: C<sub>13</sub>H<sub>20</sub>N<sub>4</sub>O<sub>2</sub>S

纯度: 98.3%

InChi Key: InChIKey=PNYYVHOTXOEBEV-UHFFFAOYSA-N

Smiles: O=C(OC1=NSN=C1N2CCCCC2)N3CCCCC3

外观: 固体粉末

作用通路: Metabolism

作用靶点: Lipase

溶解性: DMSO: 50 mg/mL (168.7 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Lalistat 2 is a selective inhibitor of lysosomal acid lipase with IC<sub>50</sub> of 152 nM and exhibits no inhibition of human pancreatic lipase or bovine milk lipoprotein lipase (up to 10 μM).