
Allopurinol Sodium

产品编号: D57045

CAS: 17795-21-0

分子式: C₅H₄N₄NaO

纯度: 98%

InChi Key: InChIKey=PTJRZVJXXNYNLN-UHFFFAOYSA-N

Smiles: [Na+].O=C1NC=NC2=C1C=NN2

外观: 固体粉末

作用通路: Immunology/Inflammation

作用靶点: ROS

溶解性: DMSO:30 mg/mL (188.6 mM);Ethanol:2 mg/mL (12.57 mM);Water:29 mg/mL (182.3 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Allopurinol Sodium is a xanthine oxidase inhibitor with an IC₅₀ of 7.82±0.12 μM.