
kb-NB77-78

产品编号: D57054

CAS: 1350622-33-1

分子式: C₁₈H₂₅NO₃Si

纯度: 99.89%

InChi Key: InChIKey=UNMWMPXUIXEQJZ-UHFFFAOYSA-N

Smiles: CC(C)(C)[Si](C)(C)Oc1cc2c(cc1)oc(=O)c1c2CCCN1

外观: 固体粉末

作用通路: Cell Cycle/Checkpoint; Metabolism

作用靶点: Serine/threonin kinase

溶解性: DMSO: 45 mg/mL (135.75 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: kb-NB77-78 is an analog of CID797718, which is a by-product of the synthesis of the parental compound, CID755673(PKD1 inhibitor).