

---

## LH846

产品编号: D57087

CAS: 639052-78-1

分子式: C<sub>16</sub>H<sub>13</sub>ClN<sub>2</sub>O<sub>2</sub>S

纯度: 98.13%

InChi Key: InChIKey=DYHAMRNAHTWYKY-UHFFFAOYSA-N

Smiles: CC1=CC2=C(C=C1Cl)N=C(NC(=O)CC1=CC=CC=C1)S2

外观: 固体粉末

作用通路: Metabolism; Stem Cells

作用靶点: Casein Kinase

溶解性: DMSO: ≥45 mg/mL

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: LH846 is a selective inhibitor of CK1δ (IC<sub>50</sub> values are 290 nM, 1.3 μM and 2.5 μM for CK1δ, ε and α); exhibits no inhibitory activity at CK2.