
PI3K-IN-46

产品编号: D57125

CAS: 304645-61-2

分子式: C₁₃H₉N₃O₃S

纯度: 97.51%

InChi Key: InChIKey=GWSGFSPSQGXVFI-UHFFFAOYSA-N

Smiles: O=C1N=C(SC1=CC=CC=CC3=NC=CC=C3C2)N

外观: 固体粉末

作用通路: PI3K/Akt/mTOR signaling

作用靶点: PI3K

溶解性: DMSO: 2.55 mg/mL (10 mM), Need ultrasonic

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: PI3K-IN-46 is a specific inhibitor of PI3Kγ.