
VX-150

产品编号: D57127

CAS: 1793080-72-4

分子式: C₂₁H₁₇F₄N₂O₇P

纯度: 95.06%

InChi Key: InChIKey=SQDQNNKQTHOQHO-UHFFFAOYSA-N

Smiles: Cc(cc(cc1)F)c1Oc(cc(C(F)(F)F)cc1)c1C(NC(C=CN1COP(O)(O)=O)=CC1=O)=O

外观: 固体粉末

作用通路: Membrane transporter/Ion channel

作用靶点: Sodium Channel

溶解性: DMSO : 125 mg/mL (242.09 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: VX-150 is a highly selective Nav1.8 inhibitor relative to the other sodium channel subtypes (>400-fold).