
Adenosine receptor antagonist 4

产品编号: D57165

CAS: 133240-06-9

分子式: C₁₀H₁₃N₃

纯度: 98%

InChi Key: InChIKey=XWWJWZJOSWSJQV-UHFFFAOYSA-N

Smiles: CCc([nH]1)nc2c1nc(C)cc2C

外观: 固体粉末

作用通路: GPCR/G Protein; Neuroscience

作用靶点: Adenosine Receptor

溶解性: DMSO: 55 mg/mL (313.86 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Adenosine receptor antagonist 4 is an adenosine receptor antagonist.