
5,7-Dichlorokynurenic acid

产品编号: D57175

CAS: 131123-76-7

分子式: C₁₀H₅Cl₂NO₃

纯度: 99.56%

InChi Key: InChIKey=BGKFPRIGXAVYNX-UHFFFAOYSA-N

Smiles: OC(c1cc(O)c2c(Cl)cc(Cl)cc2n1)=O

外观: 固体粉末

作用通路: Neuroscience

作用靶点: NMDAR

溶解性: 1eq. NaOH:100 mM; DMSO:100 mM

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: 5,7-Dichlorokynurenic acid is an NMDA receptor antagonist.