

---

## IL-15-IN-1

产品编号: D57207

CAS: 1831830-20-6

分子式: C<sub>30</sub>H<sub>36</sub>N<sub>6</sub>O<sub>4</sub>S

纯度: 99.28%

InChi Key: InChIKey=GWNFQAKCJYEJEW-UHFFFAOYSA-N

Smiles: O=C(OCC)C1=CC=CC(NC(CCCCCCSC2=NN=C(CC3=NN(C)C(C4=C3C=CC=C4)=O)N2C)=O)=C1

外观: 固体粉末

作用通路: Immunology/Inflammation

作用靶点: IL Receptor; Interleukin

溶解性: DMSO:122 mg/mL (211.55 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: IL-15-IN-1 is a selective and potent Interleukin 15 (IL-15) inhibitor, inhibiting the proliferation of IL-15-dependent cells (IC<sub>50</sub>: 0.8 μM).