
Tetraethylenepentamine Pentahydrochloride

产品编号: D57282

CAS: 4961-41-5

分子式: $C_8H_{23}N_5 \cdot 5HCl$

纯度: 98%

InChi Key: InChIKey=VGICAQBLDJAGSJ-UHFFFAOYSA-N

Smiles: Cl.Cl.Cl.Cl.Cl.NCCNCCNCCNCCN

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO: 3.72 mg/ml (10 mM, Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Tetraethylenepentamine Pentahydrochloride is a high affinity copper and zinc chelator.