
Wnt/ β -catenin agonist 2

产品编号: D57285

CAS: 943820-93-7

分子式: C₁₃H₁₂N₄O₃

纯度: 99.58%

InChi Key: InChIKey=IXPRTBBVBMTJDL-UHFFFAOYSA-N

Smiles: O=C(NCCN1N=CC=C1)C2=NOC(=C2)C=3OC=CC3

外观: 固体粉末

作用通路: Cytoskeletal Signaling; Stem Cells; PI3K/Akt/mTOR signaling

作用靶点: Wnt/ β -catenin; GSK-3

溶解性: DMSO : 45.0 mg/mL (165.3 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Wnt/ β -catenin agonist 2 activates Wnt/ β -catenin signaling and is an effective Wnt agonist.

Wnt/ β -catenin agonist 2 can be used in the study of diseases related to signal transduction and is closely related to cancer invasion and metastasis.