
Recilisib

产品编号: D57355

CAS: 334969-03-8

分子式: C₁₆H₁₃ClO₄S

纯度: 99.84%

InChi Key: InChIKey=KBEKQQJUNVQLDZ-MDZDMXLPSA-N

Smiles: C1=CC(=CC=C1CS(=O)(=O)/C=C/C2=CC=C(C=C2)C(=O)O)Cl

外观: 固体粉末

作用通路: Cytoskeletal Signaling; PI3K/Akt/mTOR signaling

作用靶点: Akt; PI3K

溶解性: DMSO:35.71 mg/mL (106.03 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Recilisib is a radioprotectant, can activate AKT, PI3K activities in cells.