
LY 78335

产品编号: D57390

CAS: 39959-66-5

分子式: C₈H₁₀Cl₃N

纯度: 99.3%

InChi Key: InChIKey=FQTXPVLCCDQRHY-UHFFFAOYSA-N

Smiles: CC(C1=C(C(=CC=C1)Cl)Cl)N.Cl

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO:<25 mg/ml(110.36 mM);H₂O:<25 mg/ml(110.36 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: LY 78335 is a phenylethanolamine-N-methyltransferase (PNMT) inhibitor.