
TNCA

产品编号: D57678

CAS: 42438-90-4

分子式: C₁₂H₁₂N₂O₂

纯度: 99.59%

InChi Key: InChIKey=FSNCEEGOMTYXKY-JTQLQIEISA-N

Smiles: O=C(O)[C@@H]1Cc2c([nH]c3ccccc23)CN1

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO < 2.16 mg/mL (10 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: TNCA is a high-affinity CaSR ligand (CaSRL) with co-agonist activity.