
Bifluranol

产品编号: D57964

CAS: 34633-34-6

分子式: C₁₇H₁₈F₂O₂

纯度: 98%

InChi Key: InChIKey=RDVXUHOSYIBGBT-ZWNOBZJWSA-N

Smiles: [C@H]([C@H](C)C1=CC(F)=C(O)C=C1)(CC)C2=CC(F)=C(O)C=C2

外观: 固体粉末

作用通路: Endocrinology/Hormones

作用靶点: Androgen Receptor

溶解性: DMSO:90 mg/mL (307.88 mM,Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Bifluranol (BX341) has anti-androgenic activity and has shown significant anti-prostatic activity in in vivo studies for the treatment of benign prostatic hyperplasia (BPH).