
TD-5471 hydrochloride

产品编号: D57973

CAS: 530084-35-6

分子式: C₃₂H₃₂ClN₃O₄

纯度: 98%

InChi Key: InChIKey=RNKVXNOAIUAJNK-JMAPEOGHSA-N

Smiles: [C@@H](CNCCC1=CC=C(NC2=CC(=C(OC)C=C2)C3=CC=CC=C3)C=C1)(O)C4=C5C(=C(O)C=C4)NC(=O)C=C5.Cl

外观: 固体粉末

作用通路: GPCR/G Protein;Neuroscience

作用靶点: Adrenergic Receptor

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: TD-5471 hydrochloride is a selective and potent long-acting human β 2-adrenergic receptor agonist for the treatment of chronic obstructive pulmonary disease (COPD).