
L-DOPA methyl ester (hydrochloride)

产品编号: D57979

CAS: 1421-65-4

分子式: C₁₀H₁₄ClNO₄

纯度: 99.92%

InChi Key: InChIKey=WFGNJLMSYIJWII-FJXQXJEOSA-N

Smiles: Cl.COC(=O)[C@@H](N)CC1=CC=C(O)C(O)=C1

外观: 固体粉末

作用通路: GPCR/G Protein; Neuroscience

作用靶点: Dopamine Receptor

溶解性: DMSO: 27.5 mg/mL (111.03 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: L-DOPA methyl ester (LDME) is an agonist of dopamine D1 receptor.