
NHS-PEG1-SS-PEG1-NHS

产品编号: D57980

CAS: 1688598-83-5

分子式: C₁₄H₁₆N₂O₁₀S₂

纯度: 99.23%

InChi Key: InChIKey=BEMZOXCMCMCTHY-UHFFFAOYSA-N

Smiles: O=C(ON1C(=O)CCC1=O)OCCSSCCOC(=O)ON2C(=O)CCC2=O

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO : 225.0 mg/mL (515.6 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: NHS-PEG1-SS-PEG1-NHS is a reversible, linker that can be used for the interaction of biomolecules with active small molecules. NHS-PEG1-SS-PEG1-NHS can be effectively used for the linkage of protein liposomes or nanoparticles.